

Nicotine

FTC

TITL

1. Distillation	Uses entire pad ✓ suspension	Uses 4 ml aliquot of extract
X. Distillate Collection	Distills to 250 ml, aliquots diluted if necessary	Distills 250 ml and dilutes to 500 ml if appropriate
3. Caustic Addition	By opening still ✓	Through separatory funnel head on still
4. No. of positions on Griffith Still	2 x 6 ✓	10
5. Still condenser	Brass tube ✓	Glass coil
6. Acid added to still	Hydrochloric	Sulfuric
7. UV spectrophotometer	Barber-Coleman (?)	Beckman DBG
2-8. Personnel	Not rotated	Rotated among assignments ✓

PM3000879465